



World-Best Ultrasound Paper,
Fitting Perfectly to World-Renowned Thermal Printers, Sony & Mitsubishi

SUPER ULSTAR Ultrasound Paper



FDA Registered
Certified to ISO 9001 / ISO 14001 / ISO 13485



DURICO DURICO C&T, Inc.

<http://www.durico.co.kr>



Why pay more for the same quality?

Compatible Papers for Mitsubishi Printer

- ▶ No need to make any adjustment on the printer for optimum printed images



Paper Grade, Size and Compatibility

| Grade | Size | DURICO | SONY | MITSUBISHI |
|-------------------------------|------------|-------------------------|-----------|------------------|
| High Glossy | 110mmX18m | ULSTAR-1100HG | UPP-110HG | K91HG/KP91HG |
| | 84mmX12.5m | ULSTAR-840HG NEW | UPP-84HG | |
| High Density | 110mmX20m | ULSTAR-1100HD | UPP-110HD | |
| | 110mmX20m | ULSTAR-1100HD mibi | | K65HM/KP65HM |
| | 210mmX25m | ULSTAR-2100HD | UPP-210HD | |
| High Density Matt Standard | 110mmX20m | ULSTAR-1100HD Matt | | |
| | 110mmX20m | ULSTAR-1100S | UPP-110S | |
| | 110mmX20m | ULSTAR-1100S mibi | | K61S/KP61S/KP61B |
| | 84mmX13.5m | ULSTAR-840S NEW | UPP-84S | |

Packing Information

- ▶ Economic Purchase for smaller Packing Unit (5 rolls in an inner box)



- ULSTAR-1100 series:** Each roll in a pouch
5 rolls in an inner box
10 inner boxes in a master carton (50 rolls)
- ULSTAR-2100HD:** Each roll in a pouch
5 rolls in an inner box
5 inner boxes in a master carton (25 rolls)
- ULSTAR- 840 series:** Each roll in a pouch
10 rolls in an inner box
10 inner boxes in a master carton (100 rolls)

DURICO Durico Coating & Technology, Inc.

■ **Seoul Office**
Banpo-daero 26-gil 21,
Secho-gu(Secho-dong, Sausohun Bldg. 2F),
Seoul 137-873, Korea
Phone: 82-2-525-8405 Fax: 82-2-525-7461
E-mail: info@durico.co.kr

■ **H/Q and Plant**
Oedap 6-gil 33, Sangju-si,
Gyeongbuk 742-320, Korea
E-mail: info@durico.co.kr

■ **U.S. Office**
Durico Imaging Inc.
18553 S.Dominguez Hills Dr.
Rancho Dominguez, CA 90220, U.S.A
Office: 714-583-7900 Direct: 714-868-1945
Fax: 714-228-1220
E-mail: sales@duricoimaging.com

■ **Europe Office**
Durico Imaging s.a.r.l.
158 Rue Diderot, 93500 Pantin, FRANCE
Office: +33(0)9 82 56 98 59,
Fax: +33(0)9 82 59 38 59,
E-mail: euro@duricoimaging.com


Direct Thermal Film for Video Printer

Durico is specialized in Self-Manufacturing Top Quality Thermal Print Media for printing Ultrasound Scanned Images.



▶ Totally different from other Inferior Compatible Papers

▶ Highly Protective UV-Cured Top Layer




Precise, High Resolution Image Reproduction

▶ First-class Colour Intensity from Highest Black Colour Density

▶ Outstanding, Bright, Sharp & Detailed Medical Image Printing available from Precise Gradation of Grey-scale

▶ High Glossy, Brilliant Photorealistic Reproduction




Long Lasting Image Durability for High Storage Stability

▶ High Humidity and Heat Resistance for Long Storage

▶ Strong Resistance to Water, Finger Prints, and Ultrasound Gel Drops


| Humidity resistance 40°C, 90%, 7Days | | Heat resistance 50°C, Dry, 7Days | |
|---|------------------|-------------------------------------|------------------|
| ULSTAR | Other Compatible | ULSTAR | Other Compatible |
| | | | |



Trouble-free Printing & Long-Term Use of Printer

▶ Easy Tearing along the Cross Direction

▶ Smooth Head Matching Performance for Continuous, Peace of Mind Printing



▶ Minimal Curling for Flat Photo Media & Easy Filing

▶ Built-in Anti-electrostatic Layer Preventing Thermal Head Damage

